

VILEX



DYNEX

ONE SYSTEM, UNLIMITED EXTERNAL FIXATION.



DYNEXX

ONE SYSTEM, UNLIMITED EXTERNAL FIXATION.



HA Coated Wires & Pins

DYNEXX is the only adult external fixation system that brings the power of Hydroxyapatite (HA) coating to both fixation options. HA coating has been clinically proven to significantly reduce loosening and pin site infection.

System Features

The DYNEXX External Fixation System is designed to address stabilization of fractures and osteotomies, rear and midfoot arthrodesis, adult and pediatric leg lengthening, and correction of bone deformity in upper and lower extremities

DYNEXX combines innovation with traditional external fixation components to deliver unlimited external fixation all in one comprehensive system.




Double-Telescopic Struts

The double-telescoping struts provide both gradual and gross adjustment. This allows the struts to be utilized in both acute or gradual correction using CORREX deformity correction software.

Reimagined Components

DYNEXX rings provide the freedom to choose between slotted or holed technology, all in one system. The concentric ring design also allows for simple sizing up or down between rings to fit patient anatomies.



ORDERING INFORMATION:

CUSTOMERSERVICE@VILEX.COM

1 (800) 521 - 5002

QSD 8.12-59 REV A

Vertrieb in der Schweiz durch:



Aeschwahrstrasse 54 info@plusortho.ch
CH-4665 Oftringen tel +41(0)62 788 2010
www.plusortho.ch fax +41(0)62 788 2011

WWW.VILEX.COM

VILEX

111 MOFFITT STREET
MCMINNVILLE, TN 37110